

HYDROJET DRY MASSAGE BED

INNOVATION IN RELAXATION



RELIEVES PAIN
BACK, JOINTS,
AND MUSCLES



POWERFUL JETS
STRONG
MASSAGE FLOW



WATER-FREE
ECO-FRIENDLY
SOLUTION



FAST RELAXATION
AS LITTLE AS
15 MINUTES



PREMIUM COMFORT
MAXIMUM
WELL-BEING

ABOUT US



Dear Clients,

We specialize in importing used equipment from closed beauty salons, wellness facilities, and medical centres in Germany, the Netherlands, and Belgium.

We then inspect and refurbish the imported equipment. Depending on its condition, this may include replacing motor carbon brushes, restoring or replacing upholstery, treating and protecting leather, or restoring the paintwork and colour of the frame.

Once the equipment has been prepared and inspected, we offer it for sale. On our scheduled delivery days, usually at weekends, we offer free delivery if the delivery location is on the route of our planned distribution run.

We provide a 6-month warranty on the equipment we sell and also offer post-warranty service.


We specialize primarily in:

- **electric massage tables,**
- **dry hydro-massage beds – Hydrojets,**
- **electric pedicure chairs,**
- **electric beauty treatment chairs,**
- **wellness loungers,**
- **wellness pedicure chairs.**

We also work with spare-parts suppliers, enabling us to provide servicing and repairs for the equipment we offer.

By prior arrangement, we can also source and import specific equipment to order according to your individual requirements.

The carriage with jets moves along the track and adapts to the body.
The person simply lies down – no other movement is required.



HOW THE SYSTEM WORKS

1 BELT / TRACK


The belt moves the carriage. The carriage is firmly attached to it and moves along the entire length.

2 CARRIAGE WITH JETS

The carriage with massage jets glides along the track and precisely adapts to the contours of the body. The jets massage through a special membrane without direct contact.


3 JETS (2x)

Two powerful water jets massage the selected areas of the body – pressure and intensity can be adjusted.



CARRIAGE MOVEMENT

The carriage moves smoothly along the track and ensures even and effective massage over the entire length of the body.



MASSAGE JETS

Water is pressurized and passed through the massage jets, which pass through the membrane and act on the body.



PRECISE MASSAGE

Targeted water jets without direct contact.



ADAPTS TO THE BODY

The carriage perfectly follows the contours of the body.



SMOOTH AND QUIET OPERATION

High-quality drive and guidance for maximum comfort.




DURABLE CONSTRUCTION

High-quality materials and construction for long service life.

What is a Hydrojet?

A Hydrojet, also known as a dry hydro-massage bed, combines the power of warm water with the comfort of a traditional massage table.

The user lies on a flexible membrane, beneath which pressurized water jets move. The body is massaged without direct contact with water – the client remains clothed and dry.



The water circulates inside the device – the massage is performed through the membrane.

Massage without changing clothes, showering, or the need for a wet area

Key Benefits

- quick and comfortable treatment, suitable even during the working day
- adjustable massage intensity, area, and duration depending on the specific model
- warmth and water pressure create an intense relaxation experience
- automated operation with minimal staff involvement
- suitable for wellness centres, fitness facilities, hotels, physiotherapy, rehabilitation, and home wellness use

Why Clients Love the Hydrojet

The combination of warmth, rhythmic movement, and water pressure creates a pleasant full-body or localized massage.

For medically certified Wellssystem models, the manufacturer states that they can be used, among other applications, to relieve muscle tension, promote blood circulation, and support relaxation.

COMFORT

The client remains dry and can stay lightly clothed.

REGENERATION

A suitable addition to recovery after sports and physical exertion.

MUSCLE

The intensity can be adjusted from a gentle to a more intensive massage.

RELAXATION

Pleasant relaxation and a reduction in feelings of stress.

Note: Specific therapeutic indications and medical applications always depend on the particular model and its certification.

Hydrojet as a Business Opportunity

A dry hydro-massage bed can operate as a stand-alone paid service or as a premium addition to a membership, wellness stay, or rehabilitation programme. A major advantage is the consistent quality of each treatment and minimal staffing requirements.

**One Device. More Clients.
Minimal Staff Involvement.**

Target Businesses

- hotels, guesthouses, and wellness centres
- fitness centres, sports clubs, and recovery facilities
- physiotherapy, rehabilitation, and orthopaedic practices
- corporate relaxation areas
- beauty and spa facilities
- private home wellness areas

For Medistream devices, the manufacturer highlights 24/7 operation, minimal staffing requirements, a small space requirement, and rapid installation and commissioning. Actual profitability depends on the purchase price, utilisation rate, price per treatment, and operating costs of the specific model.

Revenue and Return on Investment Model

The following example provides a simple guideline for operators. The figures can be adjusted according to the location and target group.



MASSAGE PRICE	10 E	15 E	20 €
Clients per day	8	10	12
Revenue / day	80 E	150 E	240 E
Revenue / 30 days	2400 E	4500 E	7200 E

Revenue is not profit. To calculate the actual return on investment, electricity, servicing, financing, marketing, space rental, VAT, and other operating costs must be deducted.

The greatest potential is achieved when the Hydrojet is used repeatedly throughout the day.

Hydrojet for Senior Care Homes

Dry water massage is an attractive addition to relaxation and recovery services in senior care homes. The client remains dry, and the treatment does not require showering. The Hydrojet can be incorporated into the facility's wellness, activation, and additional care services.

Comfortable Treatment Directly at the Facility. No wet area required • no changing clothes • minimal staff involvement

Why Is the Hydrojet an Attractive Option for Senior Care Centres?

- the treatment is performed while lying on a comfortable surface
- the client remains dry and can stay lightly clothed
- treatment duration and intensity can be individually adjusted on suitable models
- one device can be used by multiple clients throughout the day
- ideal as an additional relaxation service directly within the senior care facility
- for medical models, indications, contraindications, and the manufacturer's instructions must be observed

Important:

For seniors with specific health conditions, the suitability of the treatment must be assessed individually, following the manufacturer's instructions and the recommendations of healthcare professionals.

VZP Contributions for Massages - 2026

V každém státu přispívá zdravotní pojišťovna na zdravotní masáže a relaxaci . Některé země přispívají i na samotné pořízení zdravotnického vybavení do provozu. Je dobré se o příspěvcích jak na masáže tak i vybavení informovat dopředu, Ukážeme si jak to například funguje v České Republice

1 000 Kč

**Adults who meet the
eligibility requirements**

1 500 Kč

Policyholders with diabetes

What Is Important for a Senior Care Home


- this is not an automatic entitlement for every senior – the client must meet the eligibility requirements of the specific VZP programme
- under the screening bonus programme, massages are explicitly listed among the eligible services
- according to VZP's general terms and conditions, seniors aged 65 and over may also use alternative methods of submitting their application instead of the My VZP app
- according to VZP's general terms and conditions, seniors aged 65 and over may also use alternative methods of submitting their application instead of the My VZP app
- before implementing reimbursements, it is advisable to verify that the specific service provided and the supporting documentation comply with VZP's current requirements

Source: VZP ČR, Prevention Fund Terms and Conditions 2026. Programmes and eligibility requirements may change in subsequent years.

Economic Potential Model

Model: a senior care home with 100 residents, with a contribution of CZK 1,000 per person used entirely for massage treatments.

$$100 \text{ Seniors} \times \text{CZK } 1,000 = \text{CZK } 100,000 / \text{Year}$$



At a price of CZK 250 per massage, a contribution of CZK 1,000 corresponds to four treatments per client. For 100 clients, this model would represent 400 massages per year.

With a 20-minute treatment, this amounts to approximately 133 hours of actual massage time per year.

The model does not guarantee revenue or return on investment. It assumes 100 VZP clients who meet the eligibility requirements. And assumes that the full contribution is used for massage treatments. The calculation does not include operating costs or the purchase price of the equipment.

Why Choose a Refurbished Hydrojet?

A professional Hydrojet is a robust piece of equipment with a high purchase price when bought new.

Premium Technology at a Fraction of the Cost of a New Device

A thoroughly inspected used or refurbished unit can therefore offer a very attractive price-to-performance ratio.

Before Delivery, We Carry Out

- complete functional testing of the equipment
- inspection of the carriage system, belts, bearings, and nozzles
- inspection of the pump, seals, and water circuit
- inspection of the electronics, frequency converter, and controls
- cleaning and hygienic preparation
- test operation and verification of all available programmes

Medistream

Dry water massage designed for wellness and professional use. The user lies on a soft membrane and is massaged by two precisely controlled jets of warm water. The manufacturer states that a typical treatment lasts approximately 15–20 minutes.

Medistream 1 • Medistream 2 • Medistream 3

Some Medistream 1 and Medistream 2 Hydrojets are equipped with a cooling loop. The cooling loop is required for continuous operation, for example in spas and wellness facilities. A water supply and wastewater connection are required.

Why Choose Medistream?

- pleasant dry hydro-massage without the need for showering
- programmable nozzle movement and individual settings
- suitable for hotels, fitness centres, wellness facilities, physiotherapy, and rehabilitation
- professional construction designed for repeated use
- possibility to create an attractive self-service paid service

For used or refurbished equipment, the specific features, cooling system, controls, and technical specifications may vary depending on the year of manufacture and model version.

Medistream 1

• RELAXATION • REGENERATION • COMFORT



Advantages:

- Model with an optional cooling loop for continuous operation – can be connected to a water supply and wastewater system.
- Fully customizable massages with centimetre-level precision.
- High pressure of up to 4.5 bar.
- Metal chassis.
- We have a reliable supplier.

Disadvantages:

- No preset programmes.
- More complex to configure.
- Not often in stock.

Suitable as an entry-level option for hotels and wellness centres.

Medistream 2



Medistream can be individually customized to meet your needs.



Advantages

- **Preset programmes for quick and easy setup. Can also be individually programmed with centimetre-level precision.**
- **Optional cooling loop with water supply and wastewater connection for continuous operation.**
- **We have a reliable supplier.**
- **Commonly available in blue-grey and white versions.**
- **Pressure up to 4.5 bar.**
- **Option to purchase in instalments.**

Disadvantages

- **Plastic chassis.**

A popular choice for spas and large wellness centres. The perfect middle-ground option.

Medistream 3



ADVANTAGES

- **Luxurious illuminated design.**
- **Touchscreen display.**
- **Combined plastic and metal chassis.**
- **Very easy to operate – almost everything can be set in just 3 clicks.**
- **We have a reliable supplier.**
- **Pressure up to 4.5 bar.**

DISADVANTAGES

- **No cooling loop – no water supply or wastewater connection required.**
- **Often sold out.**

A popular luxury choice for wellness centres.

Wellsystem and Wellsystem Medical

Wellsystem Medical is described by the manufacturer as its primary Class IIa medical Hydrojet for therapeutic dry water massage. It allows individual adjustment of pressure, massage areas, massage techniques, and treatment duration.

**Professional Solution for Physiotherapy,
Rehabilitation, and Premium Wellness**

Selected Features

- Full-body or partial-body massage.
- Individual adjustment of intensity and massage areas.
- Automatic treatment without the need for manual massage by staff.
- Comfortable use while remaining clothed.
- Available in Medical and Relax versions.
- Option of mechanical or digital controls.

For the Wellsystem Medical model, the manufacturer currently specifies a maximum power consumption of 2.8 kW. The technical specifications of each individual used unit should always be verified on the manufacturer's rating plate. ([wellsystem.com](https://www.wellsystem.com))

WELLSYSTEM



**SHOPZAR – THE LARGEST SUPPLIER OF
REFURBISHED USED HYDROJETS.**

Enjoy spa-like relaxation in the comfort of your own home
at an unbeatable price.

SHOPZAR



ADVANTAGES

- Model with an optional cooling loop for continuous operation – can be connected to a water supply and wastewater system.
- Pressure up to 3.6 bar.
- Metal chassis.
- Option to purchase in instalments.
- Usually available in stock.
- Easy to repair.

Disadvantages

- Thanks to its easy-to-repair design, it has fewer massage movements and less unnecessary electronics.

Thanks to its easy-to-repair design, it has fewer massage movements and less unnecessary electronics.

Wellsystem Medical



ADVANTAGES

- Preset programmes for quick and easy setup.
- Optional cooling loop with water supply and wastewater connection for continuous operation.
- We have a reliable supplier.
- Pressure up to 4 bar.

DISADVANTAGES

- Plastic chassis.
- Often sold out.

**THE PERFECT MIDDLE-GROUND OPTION FROM
WELLSYSTEM – POWERFUL ENOUGH EVEN FOR SPA
FACILITIES**

VelusJet

Modern German Hydrojet for professional dry water massage.

VelusJet uses dual-pump technology, with higher-end versions featuring independent nozzle control. The device is designed for long-term professional operation.

VELUSJET® • Advanced • Professional

Key Features

- Dual-pump technology
- Independent nozzle control, depending on the version
- Smart Cooling with no water consumption on supported versions
- NFC / card-based control, depending on configuration
- Remote maintenance and upgrade options
- Professional: 10" touchscreen and a library of themed massage programmes
- Certified as a Class IIa medical device

Modern Hydrojet for Wellness, Fitness, Physiotherapy, and Rehabilitation

The technical specifications of each individual used unit must be verified according to its rating plate and specific equipment configuration.

VELUSJET



ADVANTAGES

- **Luxurious illuminated design.**
- **Touchscreen display.**
- **Very easy to operate – almost everything can be set in just 3 clicks.**
- **Pressure up to an INCREDIBLE 7.5 BAR from two pumps, one for each nozzle.**
- **The nozzles can be adjusted independently, with a wide selection of massage programmes available.**
- **Combined metal and plastic chassis.**
- **Optional cooling loop with water supply and wastewater connection for continuous operation.**

DISADVANTAGES

- **NO SUPPLIER AVAILABLE**
- **LIMITED NUMBER OF UNITS**
- **OFTEN SOLD OUT**
- **DIFFICULT TO REPAIR, ESPECIALLY THE ELECTRICAL SYSTEM.**

**THE MERCEDES OF HYDROJETS. IDEAL
FOR SPAS AND LUXURY WELLNESS
CENTRES.**

Installation and Operation

Requirements vary depending on the model. Models with a water cooling loop require a water supply and wastewater connection.

Because the water circulates in a closed system.

However, for each specific device, it is always necessary to verify the cooling method, electrical protection requirements, dimensions, weight, and service access.

- The buyer is responsible for arranging the unloading of the Hydrojet from the delivery vehicle and its transport to the location where it will be used.
- As our delivery vehicle is operated by a two-person team, we are unable to ensure the safe transport of the Hydrojet from the vehicle to its final installation location.
- We provide moving dollies, assistance, and our experience with handling and moving Hydrojets. However, the buyer bears full responsibility for ensuring safe and proper handling during unloading and transport to the final location.

Guaranteed 6-Month Warranty and Optional Post-Warranty Service.

Space Requirements	Check the dimensions of doorways, corridors, and the installation area.
Electricity	Verify the voltage, power consumption, and required circuit protection according to the rating plate. A B16 circuit breaker is recommended.
Cooling	Determine whether the device has air cooling or water cooling.
Service	Leave sufficient space around the device for maintenance and servicing.
Operation	Set the massage duration, price, and operating method.

Hydrojet in Your Business

Give your clients a service that combines relaxation, modern technology, and easy operation.

Hydrojet can be offered as a stand-alone paid treatment, part of a wellness package, or a benefit for fitness centre members and employees.

15-20 Minutes of Relaxation That Can Work for Your Business All Day Long.

Business Use Opportunities

- Payment per individual massage.
- Package of 5 or 10 treatments.
- Add-on to a hotel stay.
- Premium fitness membership.
- Recovery programme for athletes.
- Corporate wellness and employee benefits.

HYDROJET

Dry Water Massage for Modern Wellness

**Sales • Consulting • Technical
Inspection • Delivery • Installation**

Medistream • Wellssystem • Velus jet

Šárka Kubelková

Čs. Partyzánů 20, 537 01 Chrudim

Company ID: 05563411

Mobile: +420 720 739 640

Email: Kublkova.s@seznam.cz

SOURCES OF TECHNICAL INFORMATION

Wellssystem, Medistream, and VelusJet – official information on Hydrojet massage and individual models, including official descriptions of the technology, operating principles, and professional applications. The information has been used as a general technical reference. The specifications of each individual used unit must always be verified according to its rating plate and specific equipment configuration.

PRICE LIST

HYDROJET – Dry Massage Bed 2026

All previous price lists become invalid upon the publication of this new price list.



PRICE LIST HYDROJET

We reserve the right to change prices due to technical modifications. Delivery dates must be coordinated with us and are binding only if confirmed by us in writing.

All prices are stated excluding VAT.


Cena

4 800 E

6 000 E

10 000E

Medistream 1



Medistream 2



Medistream 3



PRICE LIST OF HYDROJETS



4 800 E



8 000 E



12 000 E